



BAXTER BARD 300XL SYRINGE INFUSION PUMPS



DIMENSIONS:	8.5 x 4.5 x 1.6 inches
WEIGHT:	1 lbs. 12 oz with batteries
ACCURARY TRAVEL:	± -3%
BATTERY LIFE:	400 hours typical minimum at low battery
BATTERIES:	Four C-size alkaline cell
FLOW PROFILE:	Continuous Flat
RATE (TRAVEL):	Xlow: 1.8 in/hr Low: 2.7 in/hr Med: 5.5 in/hr *Consult Flow Rate Chart to determine flow rate with specific syringes
OCCLUSION FORCE:	8 lbs *Consult Flow Rate Chart to determine occlusion pressure with specific syringes
STORED VOLUME ON OCCLUSION:	0.9 mL
BACK PRESSURE EFFECT ON ACCURACY:	None to occlusion pressure



844-DIAL-J2S
WWW.J2SMEDICAL.COM