

CADD LEGACY PCA AMBULATORY INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	1.6" x 3.8" 4.4"
WEIGHT	13.8 oz including 2 AA batteries, empty 100mL medication cassette reservoir
POWER SOURCE	Two AA alkaline batteries; AC adapter
MOISTURE PROTECTION	Splashproof (IPX4)
MAXIMUM INFUSION PRESSURE	40.0 psi
SYSTEM OPERATING TEMPERATURE	35°F to 104F°
SYSTEM DELIVERY ACCURACY	+ 6% nominal, mL
BOLUS ACCURACY SPECIFICATION	+6%