



## BAXTER PCA II SYRINGE INFUSION PUMPS



|                           |  |
|---------------------------|--|
| DIMENSIONS:               | 13" x 6.3" x 2.8" (330 x 160 x 71 mm)  |
| WEIGHT:                   | 4.2 pounds (1.9 kg) with batteries   |
| ACCURACY DOSE:            | ±3% average linear displacement. Rate: ±3% average linear displacement                                   |
| BATTERY LIFE:             | 30 days typical PCA and Basal/PCA Modes only. Type: four D-size alkaline cells or rechargeable batteries |
| SYRINGES:                 | Baxter® 50 mL Prefilled Syringes Monoject<br>60 cc plastic B-D 60 cc plastic                             |
| FLOW RATES:               | 0.1 to 150 mL/hr, depending on the program/cartridge   |
| OCCLUSION FORCE PRESSURE: | 9.5 ± 1.3 pounds; equivalent to 11 psig nominal  |
| INDICATORS:               | alphanumeric LCD with 8 lines x 14 characters each   |



844-DIAL-J2S  
WWW.J2SMEDICAL.COM



|                       |   |
|-----------------------|---|
| SECURITY FEATURES:    | Locking syringe cover with key lock<br>Lockable IV pole clamp Security code<br>Program cartridge contained within locking syringe cover |
| FRONT PANEL CONTROLS: | START/STOP, ENTER, HISTORY, CLEAR 0-9<br>decimal points ON and OFF switches inside locking syringe cover                                |



844-DIAL-J2S  
WWW.J2SMEDICAL.COM