



LIFECARE PCA™ PLUS II SYRINGE INFUSION PUMPS



DIMENSIONS	8" W x 13" H x 6" D
WEIGHT	15 lbs
LOCKOUT INTERVAL RANGE	5 to 100 minutes in 1-minute increments
BACKPRESSURE RANGE	10 to 20 psig
OPERATING ENVIRONMENTPOWER REQUIREMENTS	50 to 104 ° F ambient temperature
POWER REQUIREMENTS	110-120 VAC, 50/60Hz, 25 W
FUSE	0.250 amperes, time lag
BATTERY	One 8 VDC, four-cell sealed lead-acid rechargeable battery pack. When fully charged, a battery will typically provide a minimum of four hours of operation.
POWER CORD	A nominal 3-meter (10-foot) Ulrtype line power cord with integral transparent three-prong hospital-grade connector.
ELECTRICAL SAFETY	Meets UL 544 and CSA 22.2 standards
ADMINISTRATION SETS	Use only compatible Abbott PCA sets with integral Anti-Siphon valve
PCA MODE	Approx. 1 mL in 35 seconds



844-DIAL-J2S
WWW.J2SMEDICAL.COM



CONTINUOUS MODE	Variable from 0.1 x concentration (mg/hr or ug/hr) to 20.0 x concentration (mg/hr or ug/hr)
PCA + CONTINUOUS MODE	Variable from 0.1 x concentration (mg/hr or J.Lg/hr) to 20.0 x concentration (mg/hr or J.Lg/hr). + PCA dose (mg)
KVO (KEEP VEIN OPEN) RATE	0.1 mL/hr



844-DIAL-J2S
WWW.J2SMEDICAL.COM