



Hospira QCore Sapphire™ Multi-Therapy Infusion Pump Product Information

DEVICE ACCESSORIES INCLUDE:

- Pole Clamp
- AC Adapter
- cradle with splitter
- Bolus cord



Hospira Product List Number	16365
Pump Line	Sapphire™
Pump Type	Ambulatory Infusion
Dimensions Approx.	5.63" h x 3.78" w x 1.93" d
Weight Approx.	14.7 oz (without battery)
Battery Life Approx.	24 hours at 125 mL/hr
Power Requirements	100-240 V, 50-60 Hz, 0.6 A
Volume Range	1 - 9999 mL (increments of 1 mL)
Delivery Rate Accuracy	± 2.5% (subject to external conditions)
Delivery Rates	0.1 - 99.9 mL/hr with increments of 0.1 mL/hr; 100 - 999 mL/hr with increments of 1 mL/hr
System Alarms	Air In Line, Cassette Misplaced, Occlusions, Check for Occlusion, Downstream Occlusion, Upstream Occlusion, Flow Error, Insufficient Battery
Occlusion Alarm Pressure	Up to 17.4 PSI (1.2 Bar)
KVO Rate	Up to 20 mL/hr with increments of 0.1 mL/hr

All product is sold patient ready with warranty, to order call 844-DIAL-J2S
Thank you for choosing J2S Medical, we appreciate your business.