

ALARIS® MEDSYSTEM III INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	7.875" x 6" x 2.10"
WEIGHT	5.1 LBS
AIR-IN-LINE	500 ul, Neonatal 60 un
OCCCLUSION PRESSURE	15psi, Controller pressure device is 3ft H2O
OPERATING TEMPERATURE	50-104° Fahrenheit
STORAGE TEMPERATURE	Less than 95° Fahrenheit
MAXIMUM STORAGE TEMPERATURE	131° Fahrenheit



DISPLAY	Visual and audible prompts instructing user in steps to be taken in order to provide safe and appropriate therapeutic intervention
RATE RANGE	0.1-999 milliliter per hour (each channel)
VOLUME RANGE	0.1-999 milliliter (each channel)
KVO RATE RANGE	0.1-20.0 milliliter per hour
SYSTEM ACCURACY	0.1-999 ml/hr + 10% with a standard deviation of 1.96 under specified conditions
ADMINISTRATION SETS	Use only IVAC MedSystem III Administration Sets
POWER CONSUMPTION	6 Watts AC power. Use only IVAC MedSystem III AC Adapter, Model 1555 or 1550



844-DIAL-J2S
WWW.J2SMEDICAL.COM