

BAXTER FLO GUARD 6200 INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	9.50" x 6" x 11"
WEIGHT	14 LBS
LEAKAGE CURRENT	Less than 50 microamps
FUSE	0.5A, 250V, SB, 1/4in x 1-1/4in
POWER CORD	8 ft long, hospital grade plug
AC POWER REQUIREMENTS	110/120V, 60Hz
BATTERY RECHARGE	16 hours to recharged to 80%
BATTERY LIFE	5 hours on fully charges batteries



KVO RATE	5 mL/hr or programmed rate, whichever is less
VOLUME RANGE	1-9999 mL for primary and secondary
FLOW RATE RANGE	1-999 mL/hr in 1 mL increments, for primary and secondary
AIR-IN-LINE DETECTION	Air bubbles approximately 50 microliters or larger
OUTPUT PRESSURE	Initial distal occlusion alarm pressure is approximately 5 to 7 psi
ADMINISTRATION SET	Baxter standard administration sets with s suffix



844-DIAL-J2S
WWW.J2SMEDICAL.COM