

BBRAUN HORIZON NXT INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	12.4" x 4.8" x 12.2"
WEIGHT	11.9 LBS
POWER REQUIREMENTS	90-132 VAC, 30 watts, 47/63 Hz excluding accessories, 180-264 VAC, 30 Watts, 47/63 Hz excluding accessories
GROUNDING RESISTANCE	Meets UL Standard 544
LEAKING CURRENT	Meets or exceeds UL Standard 544, Section 27.2, and ANSI/AAMI Safe Current Limits
PLUG	Hospital Grade
POLE SIZE RANGE	0.75" to 1.25"
BATTERY	12 Volt sealed lead-acid battery for ambulatory use. Run time is dependent on the state of battery charge and rate of fluid delivery.



RATE RANGE	Standard: (Primary and Piggyback) 0.1 to 999.9mL/hr Micro: (Primary and Piggyback) 0.1 to 99.9 mL/hr
KEEP OPEN RATE	3 mL/hr for prescribed delivery rates of 3 to 999.9 mL/hr, or the continuation of the prescribed delivery rate for 0.1 to 2.9 mL/hr settings.
MEMORY	Permanent data retention after pump is turned off, unless battery is depleted or disconnected, or unless Data Retention Defaults have been otherwise selected in Biomed Options.
FLUID TYPES	All standard IV Fluids
HOLD TIME ALARM	Activated after three minutes without user interaction



844-DIAL-J2S
WWW.J2SMEDICAL.COM