

BBRAUN INFUSOMAT SPACE INFUSION PUMPS

ACCESSORIES INCLUDE:

- Hook and Loop tape
- Line fixation
- Metal Clip
- Cable Strap
- PCA-Key



DEVICE SPECIFICATIONS

DIMENSIONS	8.42" x 2.67" x 4.88"
WEIGHT	3.1 LBS
TIME OF BATTERY PACK	NiMH
OPERATING TIME OF BATTERY	4 hours at 100mL/hr
KVO RATE	Delivery Rate > 10ml/hr: KVO rate 3ml/hr Delivery Rate < 10 ml/hr: KVO rate 1ml/hr Delivery Rate < 1 ml/hr: KVO rate set rate



ACCURACY OF SET DELIVERY RATE	+ 5%
COMPUTER CONNECTION	USB connection in combination with BBraun interface lead CAN SP including electrical insulation



844-DIAL-J2S
WWW.J2SMEDICAL.COM