

BBRAUN OUTLOOK 100 INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	12" x 4.5" x 8.75"
WEIGHT	11 LBS
POWER REQUIREMENTS	120 VAC, 2A, 60 Hz excluding accessories
PLUG	Hospital Grade
POLE SIZE RANGE	0.75 – 1.25
BATTERY OPERATING TIME	3 hours operating time at a rate of 125 mL/hr
KVO RATE	3.0 mL/hr or selected rate, whichever is lower
VOLUME DELIVERY ACCURACY	+ 5%



FLUID TYPES	All standard IV fluids
AIR IN LINE ALARM	Greater than 100 microliter bubble
ALARMS	Air-In-Line, Battery Very Low, Close Roller Clamp, Container Empty, Downstream Occlusion, Door Open, Hold Time Exceeded, KVO, Low Battery Warning, System Error



844-DIAL-J2S
WWW.J2SMEDICAL.COM