

BBRAUN OUTLOOK 200 INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	12" x 4.5" x 8.75"
WEIGHT	11 LBS
BATTERY OPERATING TIME	3 hours
DELIVERY RATE RANGE	0.1 to 999.9 mL/hr, in 0.1 mL/hr increments
OCCCLUSION LIMIT SETTINGS	75, 100, 200, 300, 400, 500, and 750 mmHg
KVO RATE	3.0 mL/hr or selected rate, whichever is lower
POLE SIZE RANGE	0.875" - 1.3"
POWER REQUIREMENTS	120 VAC, 2A, 60 Hz excluding accessories