

BBRAUN VISTA INFUSION PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	5.5" x 9.4" x 7.8"
WEIGHT	6.8 LBS
BATTERY	7.2V, 1.2 Ah NiCd (rechargeable)
BATTERY OPERATING TIME	3.5 hours at 125 mL/hr
OCCCLUSION PRESSURE	Low: 8 PSI High: 17 PIS
KVO RATE	1mL/hr or selected rate, whichever is lower
FLOW RATE ACCURACY	+ 5% measurement time 8hr
TUBING	B. Braun Vista administration sets
ALARMS	Air-In-Line, Low Battery, Plug In Now, Bag Empty, Downstream Occlusion, Door open, Pump Stopped, KVO, Unit Malfunction