

CAREFUSION ALARIS® 8110 SYRINGE

ACCESSORIES INCLUDE:

- Disposable syringes
- Administration sets



DEVICE SPECIFICATIONS

DIMENSIONS	6.9" x 8.8" x 9"
WEIGHT	7.2 LBS
BATTERY OPERATION	Run time is a function of number of modules attached and module activity. 6 hours with one Pump Module infusing at 25 mL/h 6 hours with one Pump Module infusing at 25 mL/h and one Auto ID-Module 3 hours with four Pump Modules infusing at 25 mL/h 3 hours with four Pump Module infusing at 25 mL/h and one Auto ID-Module 4.5 hours with one active SPO2 Module 6 hours with one Syringe Module or PCA Module infusing at 5mL/h 3 hours with four Syringe Modules, or one PCA Module and three Syringe Modules, infusing at 5 mL/h 4 hours with one active EtCO ₂ Module



COMMUNICATION DATA PORT	RS-232 with an RJ45 connector
ELECTRONIC MEMORY	System configuration parameters stored in volatile memory are retained for at least six months by internal backup lithium battery. Module specific parameters are stored for 8 hours when system is off. After 8 hours of continuous off time or if a module is detached, module specific trend data and module specific operating parameters are automatically purged.
MODE OF OPERATION	Continuous
POWER REQUIREMENTS	100-240V~, 50/60 Hz, 150 VA MAX



844-DIAL-J2S
WWW.J2SMEDICAL.COM