

## CAREFUSION ALARIS® 8120 PCA SYRINGE

### ACCESSORIES INCLUDE:

- Disposable syringes
- Administration sets



### DEVICE SPECIFICATIONS

DIMENSIONS	4.5" x 15.0" x 7.5"
WEIGHT	5.5 LBS
BOLUS DOSE RANGE	Configured according to hospital best-practice guidelines
CRITICAL VOLUME	Maximum over infusion which can occur in the event of a single fault condition will not exceed 2% of nominal syringe fill volume during loading and 1% of maximum syringe travel after syringe loading.
DELIVERY UNITS	mcg, mcg/h, mg, mg/h, mL, mL/h
LOADING DOSE RANGE	Configured according to hospital best practice guidelines.



MAXIMUM DOSE RANGE	Configured according to hospital best practice guidelines.
OCCLUSION ALARM THRESHOLDS	Three settings: Low, Medium, High
OPERATING PRINCIPLE	Positive Displacement
PCA DOSE RANGE	Configured according to hospital best practice guidelines
RATE ACCURACY	+ 2% of full scale plunger travel (not including syringe variation)
SHOCK PROTECTION	Type CF, Defibrillator Proof (PCA Module) Type BF, Defibrillator Proof (dose Request Cord)



844-DIAL-J2S  
WWW.J2SMEDICAL.COM