

CODEMASTER XL DEFIBRILLATOR

ACCESSORIES INCLUDE:

- Built in recorder
- ECG
- Battery
- AC Power backup



DEVICE SPECIFICATIONS

DIMENSIONS	15.75" x 12.75" x 8.0"
WEIGHT	24 LBS
OUTPUT ENERGY	2,3,5,7,10,20,30,50,70,100,150,200,300,360 joules
WAVEFORM	Damped sinusoidal
CHARGE TIME (BATTERY OPERATION)	Less than 5 seconds to 360 joules
BATTERY TYPE	Rechargeable sealed lead-acid. 4 Ah, 12 V nominal
PADDLES	Standard paddles are anterior/anterior, adult and pediatric. Adult electrodes slide off to expose pediatric electrodes. Paddle cord is 10ft.
DISPLAY SIZE AND TYPE	5 in diagonal CRT for 4 seconds of ECG data on screen; non-fade, fixed trace. Scrolling trace is selectable.



SWEEP SPEED	25 mm/sec nominal
FREQUENCY RESPONSE	0.5 to 40 Hz
ECG OUTPUT	1 V/mV



844-DIAL-J2S
WWW.J2SMEDICAL.COM