



CADD 5400 PLUS AMBULATORY INFUSION PUMPS

ACCESSORIES INCLUDE:

- Battery (9-volt)
- 50 mL Medication Cassette Reservoir (non-sterile, for demonstration only)
- Carrying case & Carrying pouch
- Operator's Manual with warranty information



DEVICE SPECIFICATIONS

DIMENSIONS	1.1" x 3.5" x 6.4"
WEIGHT	15 oz
MEDICATION CASSETTE RESEVOIR CAPACITY	50ml or 100ml
PUMP RESOLUTION	50 microliters per interval
POWER SOURCE	9-volt alkaline battery
SYSTEM DELIVERY ACCURACY	+ 6%, normal



844-DIAL-J2S
WWW.J2SMEDICAL.COM



PUMP TIMING ACCURACY	+ 0.1%
PUMP ALARMS	Low battery, battery depleted, pump in stop mode, controller/microprocessor or motor fault, improper delivery, power up fault, low residual volume, RES VOL = 000, high delivery pressure



844-DIAL-J2S
WWW.J2SMEDICAL.COM