



## CADD 5800 PCA AMBULATORY INFUSION PUMPS

### ACCESSORIES INCLUDE:

- Alkaline battery (9-volt)
- 50 mL medication cassette reservoir (non-sterile, for demonstration only)
- Key
- Carrying case and Carrying pouch
- Operator's Manual with warranty information
- Remote DOSE cord/button with attached clip



### DEVICE SPECIFICATIONS

DIMENSIONS	1.1" x 3.5" x 6.4" including 50ml Medication Cassette reservoir
WEIGHT	15 oz including battery and empty 50ml medication Cassette Reservoir
POWER SOURCE	9 volt alkaline battery
PUMP TIMING ACCURACY	+0.1%
SYSTEM DELIVERY ACCURACY	+ 6% nominal, mL



844-DIAL-J2S  
WWW.J2SMEDICAL.COM