

KANGAROO 324 ENTERAL FEEDING PUMPS



DEVICE SPECIFICATIONS

DIMENSIONS	7.25" x 5.50" x 4.25"
WEIGHT	3.625 lbs
PUMPING MECHANISM	Rotary Peristaltic
POWER REQUIREMENTS	120 V, 60 Hz, 1 Amp
OPERATING RANGE	5-300 ml/hr in 5 ml/hr increments
ACCURACY	+10% of selected flow rate
OCCCLUSION PRESSURE	12 psi nominal
BATTERY	Rechargeable