

## PHILIPS HEARTSTART DEFIBRILLATOR

### ACCESSORIES INCLUDE:

- Battery Pack
- SMART Pads Cartridge



### DEVICE SPECIFICATIONS

DIMENSIONS	2.80" x 7.40" x 8.30"
WEIGHT	3.3 LBS
WAVEFORM	Biphasic truncated exponential waveform of 128 J, 150 J, 155 J, 157 J, 159 J, 160 J, & 158 J
CONTROL	Green On/Off button Blue i-button Orange Shock button
CHARGE TIME (BATTERY OPERATION)	Non-Rechargeable
BATTERY TYPE	HeartStart M5070A



ADULT ELECTRODES TYPE	Single use pre-attached combined ECG sensor//Defibrillation electrodes
DISPLAY	Visual and audible prompts instructing user in steps to be taken in order to provide safe and appropriate therapeutic intervention



844-DIAL-J2S  
WWW.J2SMEDICAL.COM