

# ROSS PATROL ENTERAL FEEDING PUMP



## DEVICE SPECIFICATIONS

DIMENSIONS	8.5" x 6.5" x 4.8"
WEIGHT	6.6 lbs
PUMPING MECHANISM	Rotary Peristaltic
POWER REQUIREMENTS	100 to 130 VAC, 60 Hz, 1 Amp
OPERATING RANGE	1-300 mL/hr
ACCURACY	±10% or ±0.5mL/hr, whichever is greater
OCCCLUSION PRESSURE	18 psi nominal
BATTERY	Rechargeable sealed lead acid, 6 V; the pump will fully recharge in a maximum of 4.5 hours while pump is off, or 12 hours if the pump is operating, When fully charged, 8 hours at 125 mL/hr.