

SAMARITAN PAD DEFIBRILLATOR

ACCESSORIES INCLUDE:

- Pad-Pak
- Battery
- Electrodes



DEVICE SPECIFICATIONS

DIMENSIONS	8.0" x 7.25" x 1.9"
WEIGHT	2.4 LBS
WAVEFORM	Self-Compensating Output Pulse Envelope Biphasic escalating waveform of 150 J, 150 J, & 200 J
CONTROL	On/Off and Shock
CHARGE TIME (BATTERY OPERATION)	150J in less than 8 sec



BATTERY TYPE	Lithium Manganese Dioxide 18V , 0.8 Amp Hrs
ADULT ELECTRODES TYPE	Single use pre-attached combined ECG sensor//Defibrillation electrodes
DISPLAY	Visual and audible prompts instructing user in steps to be taken in order to provide safe and appropriate therapeutic intervention



844-DIAL-J2S
WWW.J2SMEDICAL.COM