

## ZOLL M SERIES DEFIBRILLATOR

### ACCESSORIES INCLUDE:

- Service Manual
- Electrode Pads
- Replacement Batteries
- Carry case



### DEVICE SPECIFICATIONS

DIMENSIONS	6.8" x 10.3" x 8.2"
WEIGHT	11.5 LBS
OUTPUT ENERGY	2,3,5,7,10,20,30,50,70,100,150,200,300,360 joules
WAVEFORM	Dampened sinusoid
CHARGE CONTROL	Push button on apex paddle and on front panel



CHARGE TIME (BATTERY OPERATION)	4 hours or less with integral charger
BATTERY TYPE	Sealed lead acid battery pack; 5 cells, 2V/cell
PADDLES	Standard paddles are anterior/anterior, adult and pediatric. Adult electrodes slide off to expose pediatric electrodes. Paddle cord is 10 ft.
DISPLAY	5.66 inch diagonal display
SWEEP SPEED	25 mm/sec nominal



844-DIAL-J2S  
WWW.J2SMEDICAL.COM