

SIGMA 8000 INFUSION PUMP

ACCESSORIES INCLUDE:

- Tandem Carrier
- Single Pole
- Multi-pole



DEVICE SPECIFICATIONS

DIMENSIONS	9.5" x 6.5" x 6.75"
WEIGHT	10 LBS
PUMP TYPE	Volumetric Infusion Pump
BATTERY TYPE	Type: 12-volt battery pack consisting of 6 each 22.2V/2.5 amp rechargeable sealed lead acid batteries Operation: (for batteries in "as new" condition only). >4.5 hrs. to LOW BATTERY indication. >30 min. from LOW BATTERY to PLUG IN alarm. Recharge time: 18 hrs.
ANTI-FREE-FLOW SYSTEM	Set based, utilizing IV set slide clamp
MAXIMUM PUMP PRESSURE	55 PSI
OCCCLUSION PRESSURE	Adjustable: High (19 ±9 PSI), Medium (13 ±6 PSI), and Low (6 ±4 PSI)



FLOW RATE	$\pm 5\%$ from 2-999 ml/hr and $\pm 10\%$ from 0.1-1.9 ml/hr
TOTAL VOLUME	0.1-9999 m
KVO	1 ml/hr or the selected rate whichever is lower
PUMPING MECHANISM	Linear peristaltic
ALARMS AND ALERTS	air-in-line, upstream and downstream occlusion, improper IV set loading, infusion complete, low battery and plug in



844-DIAL-J2S
WWW.J2SMEDICAL.COM