

## SIGMA SPECTRUM INFUSION PUMP

### ACCESSORIES INCLUDE:

- Tandem Carrier
- 3 Pump Carriers
- Single Pole
- Multi-Pole



### DEVICE SPECIFICATIONS

DIMENSIONS	5.8" x 4.2" x 2.5"
WEIGHT	25 OZ
KVO RATE	1 mL/hr or actual rate, whichever is lower
FLOW RATE ACCURACY	0.5 – 999 mL/hr flow range in 0.1 mL/hr increments from 0.5 to 99.9 mL/hr and 1 mL/hr increments thereafter
AIR-IN-LINE ALARM	Dual beam ultrasonic detector alarms for large bubbles but allows smaller bubbles to pass
BATTERY TYPE	Lithium Ion, 1800 mA/h, 7.2V nominal
POWER CORD	9.75 feet



ENVIRONMENTAL LIMITS	Operating temperature: 60 to 90°F (15.6 to 32.2° C), 20 to 90% relative humidity non-condensing. Storage temperature: -4 to 120°F (-20 to 49°C), 10 to 90% relative humidity noncondensing
----------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



844-DIAL-J2S  
WWW.J2SMEDICAL.COM